

Patient Commission/Configuration ID: _____

Functional Unit Serial Number: _____ Date of Invoice: _____

1. Maintenance after 18 Months* (from date of invoice)

Item	☐ 16mm	☐ 20mm	Description	Potential Problem	Measure	Done	Date
4	see maintenance plan	see maintenance plan	sliding washer	wear	replacing	☐	
8	VE3771-012/26	VE3771-012/26	O-ring damper**	wear	replacing	☐	
9	GS1108-500	GS1108-500	sliding bushing**	wear	replacing	☐	
10	FE1027-01	FE1027-01	coil spring**	wear	replacing	☐	
17	VE3771-08/10	VE3771-11/10	O-ring for securing the spring unit	wear	replacing	☐	
-	SH5803-15/07	SH5805-15/18	spring unit, blue, normal, max. 15° range of motion	wear	☐ replacing	☐	
-	SH5803-15/15	SH5805-15/25	spring unit, green, medium, max. 15° range of motion	radial move of disc springs	☐ realigning disc springs with pliers	☐	
-	SH5803-10/21	SH5805-10/40	spring unit, white, strong, max. 10° range of motion			☐	
-	SH5803-10/31	SH5805-10/60	spring unit, yellow, very strong, max. 10° range of motion	noise of spring unit	☐ greasing spring unit with spray oil (article no. FT3000-15)	☐	
-	SH5803-05/63	SH5805-05/99	spring unit, red, extra strong, max. 5° range of motion			☐	
18	BR1110-L030	BR1211-L030	sliding bushing	wear	replacing	☐	

* depending on the assessment of the distributor of the custom-made product regarding the patient's usage behaviour

** is part of the functional unit

Location, Date _____

Signature of the Orthotist _____

2. Maintenance after 36 Months* (from date of invoice)

Item	☐ 16mm	☐ 20mm	Description	Potential Problem	Measure	Done	Date
1	SB9669-L0760	SB1069-L0960	bearing nut	wear	replacing	☐	
3	SF0393-02	-	feather key with pin	breakage	replacing	☐	
4	see maintenance plan	see maintenance plan	sliding washer	wear	replacing	☐	
6	SC1403-L08/1	SC1403-L10	countersunk flat head screw with hexalobular socket**	wear	replacing	☐	
7	SH0493-01	SH0493-01	plunger**	wear	replacing	☐	
8	VE3771-012/26	VE3771-012/26	O-ring damper**	wear	replacing	☐	

Patient Commission/Configuration ID: _____

Functional Unit Serial Number: _____ Date of Invoice: _____

2. Maintenance after 36 Months* (from date of invoice)

Item	☐ 16mm	☐ 20mm	Description	Potential Problem	Measure	Done	Date
9	GS1108-500	GS1108-500	sliding bushing**	wear	replacing	☐	
10	FE1027-01	FE1027-01	coil spring**	wear	replacing	☐	
11	SC1403-L08/1	SC1403-L08/1	countersunk flat head screw with hexalobular socket**	wear	replacing	☐	
14	SC1405-L12	SC1416-L14	countersunk flat head screw with hexalobular socket**	wear	replacing	☐	
15	SC1405-L12	SC1416-L14	countersunk flat head screw with hexalobular socket (axle screw)**			☐	
5-15	SH7973-AL	SH7985-AL	functional unit plug + go modularity	loss of function, wear	replacing	☐	
17	VE3771-08/10	VE3771-11/10	O-ring for securing the spring unit	wear	replacing	☐	
-	SH5803-15/07	SH5805-15/18	spring unit, blue, normal, max. 15° range of motion	wear	☐ replacing	☐	
-	SH5803-15/15	SH5805-15/25	spring unit, green, medium, max. 15° range of motion	radial move of disc springs	☐ realigning disc springs with pliers	☐	
-	SH5803-10/21	SH5805-10/40	spring unit, white, strong, max. 10° range of motion	noise of spring unit	☐ greasing coil spring with spray oil (article no. FT3000-15)	☐	
-	SH5803-10/31	SH5805-10/60	spring unit, yellow, very strong, max. 10° range of motion			☐	
-	SH5803-05/63	SH5805-05/99	spring unit, red, extra strong, max. 5° range of motion			☐	
18	BR1110-L030	BR1211-L030	sliding bushing	wear	replacing	☐	

* depending on the assessment of the distributor of the custom-made product regarding the patient's usage behaviour

** is part of the functional unit

Location, Date _____

Signature
of the Orthotist _____