

Company: ..... Customer Number: .....  
 Orthotist: ..... Date: .....  
 Patient: ..... Insurance: .....

## The Protocol for Checking the Orthosis Function is carried out:

for the current orthosis. Produced on: .....  
☐ after maintenance. Carried out on: .....  
☐ for planning a new orthosis.  
 for handing over a new orthosis. Produced on: .....

## 1. Orthosis Data

	AFO	KAFO	KO	ortho-prosthesis
Foot Piece:	short	long and partially flexible		long and rigid
Ankle Joint:	Lateral NEURO ..... no ankle joint Other: .....		Medial NEURO ..... no ankle joint Other: .....	
Knee Joint:	Lateral NEURO ..... articulated side bar KS ..... no knee joint Other: .....		Medial NEURO ..... articulated side bar KS ..... no knee joint Other: .....	

Does the orthosis comply with the configuration recommendation?      yes      partially      no

## 2. Checking the Orthosis' Alignment on the Workbench

❗ For the following sections, place the orthosis into the shoe.

2.1 The length of the orthosis' foot piece corresponds to the inner shoe length.

yes      no

2.2 The pitch of foot piece and shoe is identical.

yes      no

2.3 The toe spring is considered correctly.

yes      no


❗ Examine the orthosis' alignment laterally. If required, hold the orthosis in the position shown on the right and check the stops.

2.4 The orthosis' alignment matches the picture.

yes      no

2.5 The stops of all joints are reached.

yes      no



Basic Alignment of the Orthosis

## 3. Checking the Orthosis' Alignment on the Patient: Static

3.1 According to the configuration result, a dorsiflexion stop is recommended.

yes      no      not known

3.2 A visual stance analysis is performed.

yes

no

☐ without additional medical devices

Reason: .....

☐ with additional medical devices

.....

① For the following sections, make sure the patient is standing and wearing the orthosis and appropriate shoes.


3.3 A weight shift from one leg to the other is possible.

yes      rather yes      rather no      no

Reason: .....

.....

3.4 The patient matches the following positions the most (multiple selections possible):



free-handed  
stance possible

forward shift of the  
centre of gravity  
possible

shank vertical angle  
too small

shank vertical angle  
too wide

medical devices  
required

3.5 The stops are reached in the position ticked off at 3.4.

Ankle Joint: the dorsiflexion stop is reached.

yes      no      no dorsiflexion stop present

Knee Joint: the extension stop is not reached.

yes      no      no extension stop present



① For the following sections, always consider several step processes.  
Evaluate whether and how often the statements are true.

## 4.3 Gait Analysis: Foot

The patient touches the floor with the heel first.

always      most of the time  
sometimes      never


Reason: .....



A (passive) plantar flexion occurs.

always      most of the time  
sometimes      never

Reason: .....



The foot contact is complete.

always      most of the time  
sometimes      never

Reason: .....



A dorsiflexion occurs.

always	most of the time	sometimes
<input type="radio"/> ca. 5°	<input type="radio"/> ca. 5°	<input type="radio"/> ca. 5°
<input type="radio"/> > 5°	<input type="radio"/> > 5°	<input type="radio"/> > 5°
never		

Reason: .....



The heel lifts significantly from the ground.

always      most of the time  
sometimes      never

Reason: .....



## 4.4 Gait Analysis: Knee

The knee joint is...

flexed and the angle is... hyperextended.

☐ approx. 15°.  
☐ < 10°.      ☐ > 20°.

Reason: .....

The knee joint is...

flexed. hyperextended.

Reason: .....

The knee joint is...

flexed. hyperextended.

Reason: .....

The knee joint reaches a flexion angle of...

approx. 0°.      approx. 5°.  
< 0°      > 5°.

Reason: .....

A flexion movement occurs in the knee joint.

always      most of the time  
sometimes      never

Reason: .....

The knee joint reaches a flexion angle of...

approx. 60°.  
< 60°.      > 60°.

Reason: .....

The knee joint reaches a flexion angle of...

approx. 0°.  
< 0°.      > 0°.

Reason: .....