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Druckversion: 2020-04

Content	Page
1. Information	12
2. Safety Instructions	12
2.1 Classification of the Safety Instructions	12
2.2 All Instructions for a Safe Handling of the System Side Bars/System Anchors	12
3. Application	13
4. Scope of Delivery	13
5. Load Capacity	14
6. Tools for Mounting the System Side Bars/System Anchors	14
7. Processing	14
7.1 Fitting	15
7.2 Bending	15
7.3 Drilling	16
7.4 Processing the Surface	16
7.5 Connecting to System Case	16
8. Maintenance	16
9. Spare Parts	17
10. Accessory Parts	17
11. Disposal	18
12. CE Conformity	18
13. Legal Information	18

1. Information

These instructions for use are addressed to orthotists or qualified/trained experts and do not contain any notes about dangers which are obvious to them. To achieve maximum safety, please instruct the patient and/or care team in the use and maintenance of the product.

2. Safety Instructions

2.1 Classification of the Safety Instructions

 DANGER	Important information about a possible dangerous situation which, if not avoided, leads to death or irreversible injuries.
 WARNING	Important information about a possible dangerous situation which, if not avoided, leads to reversible injuries that need medical treatment.
 CAUTION	Important information about a possible dangerous situation which, if not avoided, leads to light injuries that do not need medical treatment.
<i>NOTICE</i>	Important information about a possible situation which, if not avoided, leads to damage of the product.

All serious incidents according to Regulation (EU) 2017/745 which are related to the product have to be reported to the manufacturer and to the competent authority of the Member State in which the orthotist or qualified/trained expert and/or the patient is established.

2.2 All Instructions for a Safe Handling of the System Side Bars/System Anchors

WARNING

Risk of Falling Due to Improper Handling

Inform the patient about the correct use of the system joint and potential dangers (e.g. breakage of the system side bar/system anchor) especially with regards to:

- moisture and water as well as
- excessive mechanical stress (e.g. due to sports, increased activity or weight gain).

WARNING

Risk of Falling Due to Permanent Higher Load

If patient data has changed (e.g. due to weight gain, growth or increased activity), recalculate the load capacity of the system joint and choose the appropriate material. For this purpose, use the Orthosis Configurator or contact Technical Support.

⚠ WARNING

Risk of Falling Due to Improper Processing

Errors in processing can lead to a breakage of the system side bar/system anchor. Bend the system side bar/system anchor as described in these instructions for use. Pay particular attention to:

- not to heat the system side bar/system anchor for bending;
- to use the drilling jig;
- that the specified bending radius is respected and
- to remove notches or residues by fine smoothing and finishing.

⚠ WARNING

Risk of Falling Due to Improper Processing with Other System Components

Errors in processing can lead to a breakage of the system side bar/system anchor or other system components. Pay particular attention to:

- fit the system side bar/system anchor correctly into the system case;
- correctly connect the system side bar/system anchor with the system case in accordance with the production technique;
- correctly connect the system anchor with the laminate of the orthosis.

NOTICE

Limitation of the Joint Function Due to Improper Processing

Errors in processing can impair the joint function. Pay particular attention to:

- grease the joint components only slightly and
- adhere to the maintenance intervals.

3. Application

System side bars and system anchors are connecting elements. A system side bar connects a system knee joint with a system ankle joint directly. The bands of a Strong Light orthosis are attached to the system side bars. A system anchor connects a system knee joint or a system ankle joint with the shell of a laminated orthosis. It is embedded into the laminate. **System side bars and system anchors** are exclusively used for orthotic fittings of the lower extremity. They must be handled by an orthotist or qualified/trained expert.

4. Scope of Delivery

Description	Quantity
system side bar/system anchor (fig. 1)	1
raised countersunk head screw with hexalobular socket (fig. 1)	2

To secure the system side bar/system anchor to the system case, use the corresponding AGOMET® adhesive that is enclosed to the system joint's scope of delivery. It can also be ordered separately (see paragraph 10).


fig. 1

5. Load Capacity

The load capacity results from the relevant patient data and can be determined by using the Orthosis Configurator. Use the system components determined by the Orthosis Configurator when producing an orthosis and mind the recommended production technique. You will find information on the production techniques in the section "Online Tutorials" on our website www.fior-gentz.com.

6. Tools for Mounting the System Side Bars/System Anchors

Tools	System Width				
	10mm	12mm	14mm	16mm	20mm
T10 hexalobular screwdriver	x	-	-	-	-
T15 hexalobular screwdriver	-	x	-	-	-
T20 hexalobular screwdriver	-	-	x	x	x

7. Processing

The system side bar/system anchor has to be customised to the patient's leg. To ensure an optimal functioning of the system side bar/system anchor, please note the following explanations on the production techniques and the processing steps.

Strong Light Technique

This production technique is characterised by its high load capacity, even though it has a lightweight construction. The system side bars can be used in orthoses with shells made of PE/PP plastics as well as in orthoses with bands made of carbon or aluminium (fig. 2).

Anchor Lamination/Prepreg Technique

In this production technique only the system anchors are embedded into the laminate. The joints' system cases remain free (fig. 3).

Joint Lamination/Prepreg Technique

In this production technique the system anchors and the joints' system cases are embedded into the laminate (fig. 4).


fig. 2


fig. 3


fig. 4

7.1 Fitting

The system side bars/system anchors are produced with oversize so that they can only be connected with the system case via press fit. For a better alignment in the system case the system side bars/system anchors are produced with a slight chamfer.

- 1 Press the system anchor to the joint's upper part by using a washer and a screw.
- 2 To do so, use a torque screwdriver and a torque of 6Nm.
- 3 Remove washer and screw again (fig. 5).


fig. 5

7.2 Bending

- Do not use a hammer to bend the system side bar.
- In order to avoid notches, use a bending iron with round edges for bending the system side bar/system anchor (fig. 6). Both bending irons with straight edges and with curved edges can easily cause breakage of the system side bar/system anchor.
- Bending is a cold working technique. Do not heat the material since the material's characteristics can change permanently.
- Do not alternate the bending direction repeatedly because this would stress the material and could cause it to break.
- To avoid fractures when bending the system anchor, make sure not to fall below the radii given in the table (fig. 7). The bending radius depends on the thickness of the material (see table).

fig. 6

Material	Calculating the Minimum Bending Radius [R*]
aluminium	$R = 11 \times \text{material thickness}$
titanium grade 2	$R = 5 \times \text{material thickness}$
titanium grade 5	$R = 10 \times \text{material thickness}$


fig. 7

* Calculation example: A system side bar made of titanium Grade 5 is 5mm thick. Multiplied by 10, the bending radius is 50mm. This value is the minimum radius.

During bending wear working clothes with long sleeves as well as gloves and goggles to avoid injuries resulting from a possible breakage of the system side bar/system anchor.

7.3 Drilling

System anchors as well as system side bars used for ankle joints are pre-drilled. For precise drilling of holes into system side bars used for knee joints as well as for system side bars available by the metre it is absolutely necessary to use the drilling jig (fig. 8). When the bores are set, fit the system side bar/system anchor as described in paragraph 7.1.


fig. 8

7.4 Processing the Surface

Before connecting the system side bar/system anchor to the system case, remove notches and other residues from the surface. Smooth and finish the surface in direction of rolling (fig. 9-10). Make sure that you do not remove too much material.


fig. 9


fig. 10

7.5 Connecting to System Case

Strong Light Technique

The system side bar must always be screwed together with the system joint or other system components and after fitting it must be adhered additionally using the enclosed AGOMET® adhesive. Secure the screws of the system side bar with a hexalobular screwdriver and with AGOMET® adhesive or LOCTITE® 243 medium strength (fig. 11-12).

fig. 11

fig. 12

Anchor Lamination/Prepreg Technique

The system anchor must always be screwed together with the system joint or other system components and after fitting it must be adhered additionally using the enclosed AGOMET® adhesive. Secure the screws of the system anchor with a hexalobular screwdriver and with AGOMET® adhesive or LOCTITE® 243 medium strength (fig. 11-12).

Joint Lamination/Prepreg Technique

The system anchor must always be sewed with the system joint or other system components (fig.13).

fig. 13

8. Maintenance

We recommend checking the system side bar for notches during a system joint's maintenance as well. In case of notches the system side bar needs to be replaced.

9. Spare Parts

Screws that are damaged during assembly or disassembly or that are worn have to be replaced with new screws of the same size.

Fastening Screw for System Side Bars and System Anchors			
Article Number	System Width	Description	Unit
SC5303-L05	10mm	raised countersunk head screw with hexalobular socket, M3 x 8	piece
SC5404-L06	12mm	raised countersunk head screw with hexalobular socket, M4 x 6	piece
SC5405-L08	14mm	raised countersunk head screw with hexalobular socket, M5 x 8	piece
SC5405-L08	16mm	raised countersunk head screw with hexalobular socket, M5 x 8	piece
SC5405-L09	20mm	raised countersunk head screw with hexalobular socket, M5 x 9	piece

10. Accessory Parts

The system joints are delivered with a small unit of AGOMET® adhesive that is sufficient for adhering the system side bar/system anchor to the system case of the system joint. If there is a need for more adhesive, you can order more.

AGOMET® Adhesive F330			
Article Number	Content	Description	Unit
KL1101	5g	set of adhesive and hardener powder, AGOMET® F330	tin
KL1100-H	30g	hardener, AGOMET® F330	tube
KL1100	800g	adhesive, AGOMET® F330	tin

11. Disposal

Dispose of the **system side bars and system anchors** properly. The product must not be disposed of with the residual waste (fig. 14). Please comply with the applicable national laws and local regulations for the proper recycling of recyclable materials.


fig. 14

For proper disposal, it is necessary to demount the **system side bars and system anchors** from the orthosis.

12. CE Conformity

We declare that our medical devices as well as our accessories for medical devices are in conformity with the requirements of Regulation (EU) 2017/745. Therefore, the FIOR & GENTZ products bear the CE marking.

13. Legal Information

With the purchase of this product, our General Terms and Conditions of Business Transactions, Sales, Delivery and Payment will apply. The warranty expires, for example, if the product is mounted several times. Please note that the product is not supposed to be combined with other components or materials than with those recommended by the FIOR & GENTZ Orthosis Configurator. Should the product be combined with products from other manufacturers, a written consent by the seller is required.

The information in these instructions for use is valid at the date of printing. The contained product information serves as a guideline. Subject to technical modifications.

Printed version: 2020-04