




Pos.	Ø	Artikelnummer	Systemgelenke	
			Systemknöchelgelenke	Systemkniegelenke
1	12 mm	GS1206-040 GS1206-045 GS1206-050 GS1206-055 GS1206-060	NEURO CLASSIC, 10 mm NEURO CLASSIC-SPRING, 10 mm NEURO CLASSIC-SWING, 10 mm NEURO VARIO-CLASSIC 2, 10 mm NEURO VARIO, 10 mm NEURO VARIO 2, 10 mm NEURO VARIO-SPRING 2, 10 mm NEURO VARIO-SWING, 10 mm NEURO SWING-CLASSIC, 10 mm NEURO SWING, 10 mm NEURO SWING 2, 10 mm	
2	12 mm	GS1207-040 GS1207-045 GS1207-050 GS1207-055 GS1207-060	NEURO SWING Carbon, 10 mm	
3	14 mm	GS1407-040 GS1407-045 GS1407-050 GS1407-055 GS1407-060	NEURO CLASSIC, 12 mm NEURO CLASSIC-SPRING, 12 mm NEURO CLASSIC-SWING, 12 mm NEURO VARIO-CLASSIC 2, 12 mm NEURO VARIO, 12 mm NEURO VARIO 2, 12 mm NEURO VARIO-SPRING 2, 12 mm NEURO VARIO-SWING, 12 mm NEURO SWING-CLASSIC, 12 mm NEURO SWING, 12 mm NEURO SWING 2, 12 mm	
4	14 mm	GS1409-040 GS1409-045 GS1409-050 GS1409-055 GS1409-060	NEURO SWING Carbon, 12 mm	
5	16 mm	GS1609-040 GS1609-045 GS1609-050* GS1609-055 GS1609-060	NEURO VARIO, 14 mm	NEURO CLASSIC zero, 12 mm NEURO VARIO zero, 12 mm NEURO CLASSIC, 12 mm NEURO VARIO, 12 mm NEURO VARIO 2, 12 mm NEURO VARIO-SWING, 12 mm NEURO ACTIVE Systemgelenkschiene, 10 und 14 mm NEURO LOCK, 14 mm
6	18 mm	GS1807-040 GS1807-045 GS1807-050 GS1807-055 GS1807-060		NEURO LOCK MAX, 12 mm NEURO FLEX MAX Sperrfunktion, 12 mm NEURO FLEX MAX Rastensperrfunktion, 12 mm
7	19 mm	GS1910-040 GS1910-045 GS1910-050 GS1910-055 GS1910-060		NEURO CLASSIC zero, 14 mm NEURO VARIO zero, 14 mm NEURO CLASSIC, 14 mm NEURO VARIO, 14 mm NEURO VARIO 2, 14 mm NEURO VARIO-SWING, 14 mm NEURO MATIC, 16 mm NEURO TRONIC, 16 mm NEURO LOCK, 16 mm
8	19 mm	GS1911-040 GS1911-045 GS1911-050 GS1911-055 GS1911-060	NEURO SWING Carbon, 14 mm	

Pos.	Ø	Artikelnummer	Systemgelenke	
			Systemknöchelgelenke	Systemkniegelenke
9	20 mm	GS2009-040 GS2009-045 GS2009-050 GS2009-055 GS2009-060	NEURO CLASSIC, 14 mm NEURO CLASSIC-SPRING, 14 mm NEURO CLASSIC-SWING, 14 mm NEURO VARIO-CLASSIC, 14 mm NEURO VARIO-CLASSIC 2, 14 mm NEURO VARIO 2, 14 mm NEURO VARIO-SPRING, 14 mm NEURO VARIO-SPRING 2, 14 mm NEURO VARIO-SWING, 14 mm NEURO SWING-CLASSIC, 14 mm NEURO SWING, 14 mm NEURO SWING 2, 14 mm	NEURO CLASSIC Carbon, 14 mm NEURO LOCK MAX, 14 mm NEURO FLEX MAX Sperrfunktion, 14 mm NEURO FLEX MAX Rastensperrfunktion, 14 mm NEURO LOCK Carbon, 14 mm
10	22 mm	GS2210-040 GS2210-045 GS2210-050* GS2210-055 GS2210-060	NEURO CLASSIC, 16 mm NEURO CLASSIC-SPRING, 16 mm NEURO CLASSIC-SWING, 16 mm NEURO VARIO-CLASSIC, 16 mm NEURO VARIO-CLASSIC 2, 16 mm NEURO VARIO, 16 mm NEURO VARIO 2, 16 mm NEURO VARIO-SPRING, 16 mm NEURO VARIO-SPRING 2, 16 mm NEURO VARIO-SWING, 16 mm NEURO SWING-CLASSIC, 16 mm NEURO SWING, 16 mm NEURO SWING 2, 16 mm	NEURO ACTIVE, 16 mm NEURO ACTIVE Systemgelenkschiene, 16 und 20 mm NEURO LOCK, 20 mm NEURO LOCK MAX, 16 mm NEURO FLEX MAX Sperrfunktion, 16 mm NEURO FLEX MAX Rastensperrfunktion, 16 mm
11	23 mm	GS2311-040 GS2311-045 GS2311-050 GS2311-055 GS2311-060		NEURO CLASSIC zero, 16 mm NEURO VARIO zero, 16 mm NEURO CLASSIC, 16 mm NEURO VARIO, 16 mm NEURO VARIO 2, 16 mm NEURO VARIO-SWING, 16 mm NEURO CLASSIC Carbon, 16 mm NEURO LOCK Carbon 16 mm
12	24 mm	GS2411-040 GS2411-045 GS2411-050 GS2411-055 GS2411-060		NEURO MATIC, 20 mm NEURO TRONIC, 20 mm NEURO HITRONIC, 20 mm NEURO LOCK MAX, 20 mm NEURO FLEX MAX Sperrfunktion, 20 mm NEURO FLEX MAX Rastensperrfunktion, 20 mm
13	24 mm	GS2413-040 GS2413-045 GS2413-050 GS2413-055 GS2413-060	NEURO CLASSIC Carbon, 16 mm NEURO SWING Carbon, 16 mm	NEURO CLASSIC Carbon, 20 mm NEURO LOCK Carbon, 20 mm
14	26 mm	GS2611-040 GS2611-045 GS2611-050 GS2611-055 GS2611-060	NEURO CLASSIC, 20 mm NEURO CLASSIC-SPRING, 20 mm NEURO CLASSIC-SWING, 20 mm NEURO VARIO-CLASSIC, 20 mm NEURO VARIO-CLASSIC 2, 20 mm NEURO VARIO, 20 mm NEURO VARIO 2, 20 mm NEURO VARIO-SPRING, 20 mm NEURO VARIO-SPRING 2, 20 mm NEURO VARIO-SWING, 20 mm NEURO SWING-CLASSIC, 20 mm NEURO SWING, 20 mm NEURO SWING 2, 20 mm NEURO HiSWING, 20 mm NEURO HiSWING R+, 20 mm	NEURO CLASSIC zero, 20 mm NEURO VARIO zero, 20 mm NEURO CLASSIC, 20 mm NEURO VARIO, 20 mm NEURO VARIO 2, 20 mm NEURO VARIO-SWING, 20 mm
15	28 mm	GS2815-040 GS2815-045 GS2815-050 GS2815-055 GS2815-060	NEURO SWING Carbon, 20 mm	
16	30 mm	GS3013-040 GS3013-045 GS3013-050 GS3013-055 GS3013-060	NEURO SWING, 24 mm NEURO SWING 2, 24 mm NEURO SWING Carbon, 24 mm NEURO HiSWING R+, 24mm	

Auf diesen Seiten sind alle Gleitscheiben in Originalgröße abgebildet. Sie können überprüfen, welche Gleitscheibe Sie besitzen, indem Sie sie auf die passende Abbildung legen. Darunter finden Sie eine Liste mit Systemgelenken, für die die Gleitscheibe eingesetzt werden kann.

* Diese Gleitscheibe wird für alle darunter aufgelisteten Systemgelenke verwendet. Für die NEURO ACTIVE Systemgelenkschiene/das NEURO ACTIVE Systemkniegelenk wird nur diese Gleitscheibe verwendet.