


# Joint Overview System Ankle Joints

Use our  
Orthosis Configurator



[www.orthosis-configurator.com](http://www.orthosis-configurator.com)



DORSIFLEXION STOP

PLANTAR FLEXION STOP

	NEURO CLASSIC	NEURO SPRING	NEURO CLASSIC-SWING	NEURO VARIO-CLASSIC	NEURO VARIO	NEURO VARIO 2	NEURO VARIO-SPRING	NEURO VARIO-SWING	NEURO SWING-CLASSIC	NEURO SWING	NEURO SWING H <sub>2</sub> O
static	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗
dynamic due to variable spring force	✗	✗	✗	✗	✗	✗	✗	✗	✓	✓	✓
dorsiflexion angle adjustable on the system joint	✗	✗	✗	✓	✓	✓	✓	✓	✓	✓	✓
dorsiflexion angle adjustable by filling the system stirrup	✓	✓	✓	✓	✓	✗	✓	✗	✗	✗	✗
static	✓	✗	✓	✗	✓	✓	✗	✓	✓	✓	✗
dynamic due to adjustable spring force	✗	✗	✓	✗	✗	✗	✗	✓	✗	✓	✓
plantar flexion angle adjustable on the system joint	✗	✗	✓	✗	✓	✓	✗	✓	✗	✓	✓
plantar flexion angle adjustable by filling the system stirrup	✓	✗	✗	✗	✓	✗	✗	✗	✓	✗	✗
water-resistant 	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓