


Joint Overview System Knee Joints

Use our
Orthosis Configurator



www.orthosis-configurator.com



	NEURO CLASSIC zero	NEURO VARIO zero	NEURO CLASSIC	NEURO VARIO	NEURO ACTIVE	NEURO MATIC	NEURO TRONIC	NEURO LOCK	NEURO LOCK MAX	NEURO FLEX MAX
CONTROL MECHANISM	free moving without posterior offset	✓	✓	✗	✗	✗	✗	✗	✗	✗
	free moving with posterior offset	✗	✗	✓	✓	✓	✓	✓	✓	✓
	automatic stance phase control (mechanical)	✗	✗	✗	✗	✓	✗	✗	✗	✗
	automatic stance phase control (electromechanical)	✗	✗	✗	✗	✗	✓	✗	✗	✗
EXTENSION STOP	permanently locked	✗	✗	✗	✗	✓	✓	✓	✓	✓
	damped	✓	✓	✓	✓	✗	✓	✓	✗	✗
	adaptable by exchanging the extension stop	✗	✓	✗	✓	✓	✗	✓	✓	✓
FLEXION STOP	adaptable by exchanging the joint's upper part	✓	✗	✓	✗	✓	✓	✗	✗	✗
	adaptable by inserting a flexion stop	✗	✗	✗	✗	✓	✗	✗	✗	✓
LOCK FUNCTION	adaptable by exchanging the extension stop and adjusting the joint's lower part by filing	✗	✗	✗	✗	✗	✗	✓	✓	✗
	adaptable by exchanging the extension stop and adjusting the flexion stop disc by filing	✗	✗	✗	✗	✗	✗	✗	✗	✓
	step lock function	✗	✗	✗	✗	✗	✗	✗	✗	✓