

CROSSROADS

Orthosis shoe. Different.

with Elastic Outsole Made of EVA Material



Size: 36-48

Material: High Quality Leather
Certified cow leather used for the upper material and calf leather used for the inner lining ensure the highest wearing comfort.

Shoe Designs: REGULAR left and right
SLIM left and right

REGULAR left/SLIM right
SLIM left/REGULAR right

CROSSROADS

Orthosis shoe. Different.

with Elastic vibram® Flat Sole



Shoe Measurements

toe spring	18mm
pitch	15mm
heel height	25mm

PARKSTREET

Orthosis shoe. Different.

with Heelless Rubber Outsole



Size: 32-48

Material: Sturdy Fabric
The sturdy fabric made of densely woven cotton is durable, easy to clean and guarantees a carefree use in everyday life.

All orthosis shoes are sold in pairs for each shoe size.

URBANSTREET

Orthosis shoe. Different.

with Heelless Rubber Outsole



Shoe Measurements

toe spring	10mm
pitch	0mm
heel height	20mm