

# NEURO SWING FIT TEST-AFO



The NEURO SWING FIT TEST-AFO with premounted NEURO SWING Carbon system ankle joint is a preproduced orthosis made of carbon fibre material. It is used as a test orthosis prior to a treatment with a custom-made AFO equipped with a system ankle joint with dynamic stops (e.g. NEURO SWING).

## The NEURO SWING FIT TEST-AFO:

- is available as a set for the right and left leg.
- is equipped with a premounted NEURO SWING Carbon system ankle joint.
- is delivered with pre-cut padding material, with self-adhesive hook dots, for the lower leg shell and one cloth bag per orthosis.
- has a hook-and-loop fastener at calf height which provides high stability.


Spring units (see back page) have to be ordered separately.

Shoe Size*	Orthosis Size	Article Number for System Width			
		12mm	14mm	16mm	20mm
29-31	XXS	AC5201-XXS/LR	AC5202-XXS/LR	-	-
32-34	XS	-	AC5202-XS/LR	-	-
35-37	S	-	-	AC5203-S/LR	-
38-40	M	-	-	AC5203-M/LR	-
41-43	L	-	-	-	AC5205-L/LR
44-46	XL	-	-	-	AC5205-XL/LR

\* Determine the shoe size and select the corresponding orthosis size. For patients with a heavy body weight or of sturdy build, it may be necessary to select a larger orthosis size.



The NEURO SWING FIT TEST-AFO should be worn with a knee stocking or a knee high compression stocking.




## Spring Unit Case



### Spring Unit Case

Item	Description	Article Number
1	1 x spring unit set in a case with tools	SF5800-C
w/o fig.	1 x spring unit case for filling by yourself with tools (without spring units)	SF5800-C/0
1a	2 x of each blue, green, white, yellow and red spring unit for NEURO SWING Carbon 12mm	see table spring units
1b	2 x of each blue, green, white, yellow and red spring unit for NEURO SWING Carbon 14mm	see table spring units
1c	2 x of each blue, green, white, yellow and red spring unit for NEURO SWING Carbon 16mm	see table spring units
1d	2 x of each blue, green, white, yellow and red spring unit for NEURO SWING Carbon 20mm	see table spring units
1e	1 x hexagonal screwdriver with spherical head, 4 x 100mm	WZ5112-4010
1f	1 x hexagonal screwdriver with spherical head, 5 x 100mm	WZ5112-5010



### Spring Units (You Need Two of the Following Spring Units for One NEURO SWING Carbon System Ankle Joint):

Unit	Description	Article Number for System Width			
		12mm	14mm	16mm	20mm
piece	spring unit, blue, normal, max. 15° range of motion	SF5801-C/15/03	SF5802-C/15/05	SF5803-C/15/07	SF5805-C/15/18
piece	spring unit, green, medium, max. 15° range of motion	SF5801-C/15/06	SF5802-C/15/11	SF5803-C/15/15	SF5805-C/15/25
piece	spring unit, white, strong, max. 10° range of motion	SF5801-C/10/12	SF5802-C/09/16	SF5803-C/10/21	SF5805-C/10/40
piece	spring unit, yellow, very strong, max. 10° range of motion	SF5801-C/10/19	SF5802-C/10/29	SF5803-C/10/31	SF5805-C/10/60
piece	spring unit, red, extra strong, max. 5° range of motion	SF5801-C/05/33	SF5802-C/05/53	SF5803-C/05/63	SF5805-C/05/99