

# Joint Overview System Ankle Joints

Use our  
Orthosis Configurator



www.orthosis-configurator.com



	NEURO CLASSIC with free moving system stirrup	NEURO CLASSIC with static system stirrup	NEURO CLASSIC-SPRING	NEURO SPRING	NEURO CLASSIC-SWING	NEURO VARIO-CLASSIC	NEURO VARIO-CLASSIC 2	NEURO VARIO	NEURO VARIO 2	NEURO VARIO-SPRING	NEURO VARIO-SPRING 2	NEURO VARIO-SWING	NEURO SWING-CLASSIC	NEURO SWING	NEURO SWING 2	NEURO SWING H2O
converting option within the modular joint group	✗	✗	✗	✗	✗	✗	✓	✗	✓	✗	✓	✓	✗	✓	✓	✗
none	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
static	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗
dynamic due to adjustable spring force	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✓	✓	✓
dorsiflexion angle adjustable on the system joint	✗	✗	✗	✗	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
dorsiflexion angle adjustable by filing the system stirrup	✗	✓	✓	✓	✓	✓	✗	✓	✗	✓	✗	✗	✗	✗	✗	✗
none	✓	✗	✓	✗	✗	✓	✓	✗	✗	✗	✓	✗	✗	✗	✗	✗
static	✗	✓	✗	✗	✓	✗	✗	✓	✓	✗	✗	✓	✓	✓	✓	✗
dynamic due to adjustable spring force	✗	✗	✗	✗	✓	✗	✗	✗	✗	✗	✗	✓	✗	✓	✓	✓
plantar flexion angle adjustable on the system joint	✗	✗	✗	✗	✓	✗	✗	✓	✓	✗	✗	✓	✗	✓	✓	✓
plantar flexion angle adjustable by filing the system stirrup	✗	✓	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗	✓	✗	✗	✗
integrated noise reduction	✗	✗	✗	✗	✗	✗	✓	✗	✓	✗	✓	✓	✗	✗	✓	✗
water-resistant	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓

DORSIFLEXION STOP

PLANTAR FLEXION STOP