FIOR & GENTZ

Therapeutic Shoes
Sold From Inventories

We're clearing our inventory!

Selected therapeutic shoes now 50% off!*

Take advantage of this limited-time offer and secure selected therapeutic shoes at great prices. The articles listed below will be permanently discontinued. Take advantage of the 50% discount now – while stocks last!

Please note: Purchased shoes cannot be exchanged or returned.

* instead of your personal discount





HEIDELBERG

Bandage Shoes
Flexible sole with pitch, 30mm thick

Article number: TS4552

HEIDELBERG has a variable volume. The partially flexible upper fabric of the shoe is elastic and adapts to different space demands (e.g. due to varying or thick bandages).

Shoe Size: 35–50



DRESDEN WIDE

Bandage Shoes
Rigid sole with extra-strong rocker bottom and pitch, 30mm thick

Article number: TS5411

DRESDEN, especially the WIDE version, provides much space for specific demands. The firm outer fabric can be separated from the soft inner fabric in order to avoid possible pressure points and prevent the manifestation of oedema. Furthermore, DRESDEN is an appropriate TWINSHOE and thus height compensation for the DIABETES therapeutic shoe WETZLAR.

Shoe Size: 35–48



KASSEL STANDARD

Bandage Shoes
Rigid sole with extra-strong rocker bottom and pitch, 30mm thick

Article number: TS5421

KASSEL, especially the WIDE version, provides much space for specific demands. The firm outer fabric can be separated from the soft inner fabric in order to avoid possible pressure points and prevent the manifestation of oedema. If applicable, a protective forefoot cap can be attached to the toe area that keeps the foot dry and warm.

Shoe Size: 35–48

NÜRNBERG STANDARD

Bandage Shoes
Rigid sole with extra-strong rocker bottom and pitch, 30mm thick

Article number: TS4311



NÜRNBERG WIDE

Bandage Shoes
Rigid sole with extra-strong rocker bottom and pitch, 30mm thick

Article number: TS4411

NÜRNBERG, especially the WIDE version, provides much space for specific demands. Besides, the firm outer fabric can be separated from the soft inner fabric in order to avoid possible pressure points and prevent the manifestation of oedema.

Shoe Size: 35–48



HAMBURG

Bandage Shoes
Rigid sole with extra-strong rocker bottom and straight top, 30mm thick

Article number: TS6121

HAMBURG has an open forefoot area with a low framing brim. If necessary, a protective forefoot cap can be attached to keep the foot dry and warm. With the convenient hook and loop fastener in the heel area, the shaft width can be adapted. Furthermore, HAMBURG is an appropriate TWINSHOE and thus height compensation for the relief shoe KÖLN.

Shoe Size: XS–XXL



BERLIN

Forefoot Relief Shoes
10° wedge sole (50:50) without forefoot support, 40mm thick

Article number: FO1002/FO1012

BERLIN has a shortened wedge sole for complete pressure relief in the forefoot area. BERLIN with protective stirrup offers additional protection against injuries (e.g. impacts). The protective stirrup is both adjustable in length and thermoformable and can thus be individually adjusted.
If necessary, a protective forefoot cap can be attached to the insole to keep the foot dry and warm. BERLIN can also be equipped with a protective stirrup (FO1000-Z5,-Z6,-Z7,-Z8) subsequently.

Shoe Size: S–XL



FRANKFURT

Forefoot Relief Shoes
Rigid sole with extra-strong rocker bottom and straight top, 40mm thick

Article number: FO7101

FRANKFURT has an open toe area with a strap. If necessary, a protective forefoot cap can be attached to keep the foot dry and warm. With the convenient hook and loop fastener in the heel area, the shaft width can be adapted.

Shoe Size: XS–XXL



KOPENHAGEN

Forefoot Relief Shoes
Rigid sole with extra-strong rocker bottom and straight top, 40mm thick

Article number: FO7100

KOPENHAGEN has an open toe area with a strap. If necessary, a protective forefoot cap can be attached to keep the foot dry and warm. With the convenient hook and loop fastener in the heel area, the shaft width can be adapted.

Shoe Size: XS–XXL



STENDAL

Forefoot Relief Shoes
10° wedge sole (50:50) with forefoot support, 30mm thick

Article number: FO1160

STENDAL has an open toe area with a strap. If necessary, a protective forefoot cap can be attached to keep the foot dry and warm. With the convenient hook and loop fastener in the heel area, the shaft width can be adapted.

Shoe Size: XS–XXL

POST-OP



STOCKHOLM

Forefoot Relief Shoes
8° wedge sole (60:40) with forefoot support, 30mm thick

Article number: FO1120

STOCKHOLM has an open forefoot area with a low framing brim. If necessary, a protective forefoot cap can be attached to keep the foot dry and warm. With the convenient hook and loop fastener in the heel area, the shaft width can be adapted.

Shoe Size: XS–XXL



MÜNCHEN

Heel Relief Shoes
5° wedge sole (20:80) without heel support, 30mm thick

Article number: FO2121

MÜNCHEN provides an open forefoot area with a low framing brim and an open heel area. If necessary, a protective forefoot cap and/or a protective heel cap can be attached.

Shoe Size: S–XL



WUPPERTAL

Heel Relief Shoes
5° wedge sole (20:80) without heel support, 30mm thick

Article number: FO2000

WUPPERTAL has an open toe and heel area and can be adapted to different needs by using the different hook and loop fasteners.

Shoe Size: S–XL

You can find more information on all our therapeutic shoes in our **Product Catalogue Therapeutic Shoes**, available in the download area of the FIOR & GENTZ website.

PR0296-GB-2026-10

FIOR & GENTZ
Dorette-von-Stern-Straße 5
21337 Lüneburg (Germany)

+49 4131 24445-0
+49 921 95659554

info@fior-gentz.de
www.fior-gentz.com